

University of Wollongong



Out-Of-Field Dose Equivalents Delivered by Passively Scattered Therapeutic Proton Beams

Anatoly Rosenfeld, PhD
Professor and Director
Centre for Medical Radiation Physics

Andrew Wroe^{1,2}, Ben Clasie³, Hanne Kooy³, Jay Flanz³, Harald Paganetti³
Reinhard Schulte²



¹ Centre for Medical Radiation Physics, University of Wollongong, Australia

² Department of Radiation Medicine, Loma Linda University Medical Centre, USA

³ Department of Radiation Oncology, Massachusetts General Hospital, USA



Background

- Protons are a useful tool for cancer treatment
- The Proton depth-dose distribution IMRT dose



Neutron Production

Neutrons with
 $E < 10 \text{ MeV}$ →



← Neutrons with
 $E > 10 \text{ MeV}$

[1] M. Moyers et al., "Leakage and scatter radiation from a double scattering based proton beamline", Med. Phys., 35, 1, 2008, pp. 128-144

Dose Equivalent Measurements

- Physical measurements in such cases are scarce
- Detection of a range of particles is required
- High spatial resolution is required for measurements within phantoms in close proximity to the field edge
- Consideration of mixed particle quality factor for the determination of dose equivalent is also difficult
- Clinically relevant measurements are required
- Measurements are necessary to validate Monte Carlo


Equipments

- Bonner sphere



LU DLUM neutron ball cart model 42-5
(diameter: bare, 2", 3", 5", 8", 12")

detector:
LiI(Eu) scintillator, 4mm x 4mm




Courtesy of Dr Matsufuji, HIMAC

- rem counter



ALOKA neutron survey meter TPS-451C


detector:
 ${}^3\text{He}$ proportional counter




- Sv (up to 15 MeV)
- energy spectra not available

Experiment (HIMAC, BIO beam line)

Courtesy of Dr Matsufuji, HIMAC



- incident beam (~160mm range in water)
 - ^{12}C -290MeV/n, ^1H -160MeV
 - 100mm dia., SOBP 60mm
- target
 - water (390mm)
- position
 - 0, 15, 30, 45, 60, 90 deg
 - 173 cm from target



veto slipped


veto fired

Research Aim

To conduct microdosimetry measurements of external field dose equivalents for passively delivered clinical treatment fields in proton therapy



What is Microdosimetry

- When measuring cellular effects, it makes sense to use detectors the same size as a cell - **MICRODOSIMETRY**
- Size = $10\mu\text{m}$ = 1/2500 inches
- Each event interacting with the volume is measured and recorded
- Suitable for measuring neutron dose
- “Dose Equivalent” can be determined




Solid State Microdosimetry

- SOI microdosimeters developed by Centre for Medical Radiation Physics at the University of Wollongong
- Small size allows for accurate measurements near the treatment field edge
- Tested extensively in Proton and Heavy Ion Therapy




Silicon-On-Insulator (SOI) Microdosimeter

- Provides true microscopic Sensitive Volumes
- Volume Area $30 \times 30 \mu\text{m}^2$ or $100 \times 100 \mu\text{m}^2$
- Volume thickness either 2, 5 or $10 \mu\text{m}$
- Array size 50 – 5000 independent cells




Experimental Setup



Dose Equivalent Determination

- Dose is determined from the $f(E)/E$ Spectra in the following method:




$$D_{Si} = \frac{\int_0^{\infty} f(E).E.dE}{M} = \frac{\int_0^{\infty} f(E).E.dE}{\rho_{Si} V n_{cells}}$$

$$\therefore D_{Si} = D_{TE} \frac{S_{Si}}{S_{TE}}$$

Dose Equivalent Determination

- The lineal Energy Spectra is determined by dividing the energy by the mean chord length $\langle l \rangle$:




$$\langle l \rangle = \frac{4V}{S\xi}$$

$$y = \frac{E}{\langle l \rangle}$$

- Where $\langle l \rangle$ in this case is $19.05\mu\text{m}$ and $\xi=0.63$ is the TE conversion factor

Dose Equivalent Determination

- A normalised dose weighted lineal energy spectra can be obtained using the following relationship:




$$d(y) = \frac{yf(y)}{y_f}$$

$$y_f = \int_0^{\infty} yf(y)dy$$

$$\int_0^{\infty} d(y)dy = 1$$

Dose Equivalent Determination

- The final step in determining dose equivalent (Sv) is to convolve the dy spectra with a quality spectra $Q(y)$.



$$H = \int_0^{\infty} Q(y) D(y) dy$$


$$\therefore H = \int_0^{\infty} Q(y) D_{TE} d(y) dy$$

Dose Equivalent Determination

- Advantages include:
 - Small Size (precision in placement)
 - Established $Q(y)$, ICRU 40
 - $Q(y)$ changes with lineal energy
 - Wide range of lineal energies (1.0-1000keV/ μm)
- Possible errors include:
 - Noise threshold removing some signal (limit 0.8 keV/ μm)
 - $Q(y)$ determined in-vitro (artificial case?) not in-vivo (real case?)
 - Assumptions in average stopping power conversion
 - Secondaries produced within device (I.Cornelius et al, **CMRP**)

Clinical Treatment Fields

- Wide range of clinical treatment fields chosen
- All parameters were reproduced from clinical settings




Clinical Disease	Beam Range in Water (cm)	Beam Modulation (cm)	Maximum Field Dimensions (cm)	Pre-collimation field size (cm ²)
Prostate	28.8	10.4	6.9, 7.7	184
Medulloblastoma (brain)	15.8	15.9	15, 16.7	547
Medulloblastoma (spine)	9.6	7.4	16.2, 6.5	366
Stereotactic (brain)	18	0	2.0 cm diameter	-
Ocular Melanoma	2.7	2.5	1.1cm diameter	-

Results

- mSv d at late field e
- 8.3 mS
- Sharpe phanto
- Q_{Avg} in
- Q_{Avg} c beam
- H simil patient surfac
- H high axis





1cm WED
 3cm WED
 4cm WED



Results

Scanning parallel to the beam at 5cm offset





- $5.5 \text{ mSv/Gy} < H_{\text{aperture}} < 0.9 \text{ mSv/Gy}$
- $7.1 \text{ mSv/Gy} < H_{\text{block}} < 0.5 \text{ mSv/Gy}$
- H_{aperture} has a different dependence on depth than H_{block}
- Scattered primary protons affects H and the determination of Q up to 22.3 cm depth (see our poster [B.Clasie et al.](#))
- Downstream of the Bragg peak, difference in H is due to n generated in the phantom



Results: medulloblastoma

Measurements lateral to the primary field

- $10 \text{ mSv/Gy} < H_{\text{cranial}} < 0.6 \text{ mSv/Gy}$ at phantom surface
- Cranial medulloblastoma results similar to those for prostate at phantom surface
- In comparison with the prostate cancer field: \downarrow Proton Energy \uparrow Field Size
- Reducing energy in the spinal case produces less dose equivalent lateral to the primary field
- Q_{Avg} is larger for the spinal case, which indicates that there is less contribution from scattered primary protons




Results: downstream of the Bragg peak



- H and Q_{Avg} on central axis is similar for both treatment configurations
- $2.3 \text{ mSv/Gy} < H < 0.6 \text{ mSv/Gy}$
- Q_{Avg} constant at 5



Results: ocular melanoma





- 27 mm range and 25 mm SOBP
- 1.1 cm diameter field

Results: ocular melanoma

- Greater than order of magnitude difference in H between fields due to lower proton energy
- Q_{Avg} is similar at close to the field and increasing with lateral distance



Results: Stereotactic beam



2 cm field diameter in MLC and brass aperture, parallel to the beam at 5cm from the field edge




- H higher for the mMLC case by as much as 50%
- Q_{Avg} mMLC is generally lower indicating leakage protons

3D SOI silicon microdosimetry: new design, CMRP

3D silicon cell array: fabricated at SNFF , UNSW, Australia, Prof A.Dzurak

Detector Array Design_3(Revised) (All dimensions are in microns)

Number of cells: 30x30
Detector Array dimension: 1439umx1439um
Bond Pad dimension: 100um x 100um
Diameter of each cell: 25.8 um
Spacing between cells: 22 um




3D SOI silicon microdosimetry: new design, CMRP

Response of new 3D SOI microdosimeter on $1\text{ }\mu\text{m}$ diameter 3 MeV alpha particles scanning microbeam (ANSTO Dr M.Reinhard).

Each cell has sensitive volume with a diameter of $6\text{ }\mu\text{m}$ and pitch $20\text{ }\mu\text{m}$

Collaboration with ANSTO heavy ions micro beam probe, measurements were done by PhD student Ms Amy Ziebell, CMRP, Uni of Wollongong




Overall Conclusion

- This work has provided an assessment of dose equivalent for clinical proton therapy configurations
- The results highlight that neutrons are present but are not in great abundance
- External field dose equivalent dose depend on the treatment situation and less with reduced proton energy
- However, the external field dose equivalent is not largely field size dependant as initially hypothesized
- Results obtained are analogous to leakage from **MLC's in IMRT or less.**

Overall Conclusions

- Leakage radiation for 6MV X-ray measured along patient central axis (depth=10 cm)
- X-Ray data normalised to $D_{\max} = 1.6$ cm
- Proton data for cranial medulloblastoma field ($15 \times 17 \text{ cm}^2$) at a WED 8.4 cm
- Proton data normalised to dose at isocenter (i.e. central axis of patient)
- Proton results still consider Q factor
- Proton results a factor of 2-3 lower than that for a $10 \times 10 \text{ cm}^2$ conventional field
- Proton results a comparable to a $10 \times 10 \text{ cm}^2$ conventional field at lateral displacement of 30 cm or greater



[2] E. Klein et al., "Peripheral doses from pediatric IMRT", Med. Phys., 33, 7, 2006, pp. 2525-2531

Are out-of-field doses less of an issue for protons?

Acknowledgements

- Physics Staff at Massachusetts General Hospital, Burn Proton Therapy Centre
- Physics Staff at the Australian Nuclear Science and Technology Organization (ANSTO)