

Patient immobilization and position verification strategies:

MDACC Experience

X. Ronald Zhu, Ph.D.
The University of Texas
M.D. Anderson Cancer Center

PTCOG 47
Jacksonville, FL
May 22-24, 2008

THE UNIVERSITY OF TEXAS
MD ANDERSON
CANCER CENTER
Making Cancer History™

Acknowledgements

Therapists

- Chuck Merrifield and the team

Physicists

- Richard Amos, MS
- Narayan Sahoo, PhD
- Lei Dong, PhD
- Richard Wu, MS
- James Yang, PhD
- Michael Gillin, PhD
- Al Smith, PhD
- Many Others

Dosimetrists

- Beverly Riley and the team

Radiation Oncologists

- Andrew Lee, MD
- Shiao Woo, MD
- James Cox, MD
- Many Others

Introduction

- Each beamline at MDACC/PTCH has three x-ray systems.
- Daily image guided treatment delivery using 2 of the 3 x-ray systems.
- Three Hitachi couch extensions: short, medium and long.
- One Head and neck extension.

Hitachi gantry with orthogonal X-ray system


Cage Source

Nozzle Detector

Nozzle Source

Cage Detector

THE UNIVERSITY OF TEXAS
MD ANDERSON
CANCER CENTER
PROTON THERAPY CENTER




Hitachi Couch Extensions




Short
84 cm

Medium
115 cm

Long
136 cm




Head & Neck Extension (HNE)



- No beam going through the edge of the couch extension

BOS headframe




WET of Couch Extension

- Everything in the beam, we need to know what is the WET and If TPS can predict it.
- Water equivalent thickness (WET) can be determined through measurements of depth dose curves with and without the extension:
 - WET - Hitachi couch extensions = 1.1 cm
 - WET (HNE) = 0.7 cm

Narayan Sahoo, PhD


WET by the Planning System

- Ranges calculated for a PA beam with and without the extension:
 - With - Use CT values of the image
 - Without - Over-writing the HNE as air (HU = -1000).



Calculated WET = 0.7 cm.

IGRT Work Flow



Imaging Room

X-ray
Controls


PIAS

X-ray



Physician Work Room

Mosaiq Screen

G2 PIAS Monitor




Prostate



Prostate immobilization:


- Rectal balloon – 65 cc Water
- Leg immobilization
- Full bladder – checked by bladder scanner

AP
DRR



AP
X-ray

Lateral
DRR



Lateral
X-Ray


Localization Results - Prostate

- Translational only
- Manual alignment
- Retrospective study
 - 3 mm geometric margin is sufficient*

* Choi *et al.* PTCOG-47, 2008

Beyond Prostate


- Two issues for posterior beam:
 - CT couch is not the same as treatment couch
 - Imaging artifacts of thin high density plate






The measured WET is 1.5 cm WET, but CT number as low as -400 HU.

Digital couch

- Digital couch is necessary:
 - Replacing CT couch with treatment couch
 - Remove imaging artifacts




Dosimetric Impact



Water phantom with and without
Digital couch




Heng Li, PhD

Craniospinal Irradiation




- Styrofoam to raise the patient from table.
- Mouth piece attached to mask to reduce rotation.


Craniospinal Irradiation



AP
DRR




AP
X-ray



Lateral
DRR



Lateral
X-Ray





Thoracic Patients



- T-bar + wing board with vacuum-bag.
- Limit the vacuum-bag above the treatment fields – skin folding*

*Richard Amos



Summary

- Immobilization is one of most important aspects in proton therapy:
 - Avoid use high density thin plate and sharp edges
 - Range uncertainty introduced by immobilization
- Localization – 2 orthogonal x-ray images:
 - Translational only
 - Manual alignment
 - Improvement in imaging alignment method